

TruSystem™ 7000 U14 OR Table

Universal Operating Table for Greater Patient Access



TruSystem™ 7000 U14 Mobile Operating Table

Make your operating room more flexible

- All the advantages of the TruSystem™ 7000 OR Table but with a shortened back section, allowing greater access for the anesthesia provider regardless of the surgery.
- The modular design along with a comprehensive line of accessories and components provides the ability to configure the table for a wide range of applications and offers surgeons greater visibility in applications such as gynecology, urology and orthopedics.
- Simplified configuration and setup with improved coupling point system.
- Radiolucent carbon fiber components are also available for full patient intraoperative imaging.

Components & Accessories for Various Applications



Superior access during urology and gynecology procedures with stirrups and URO cut-out.



Simplified patient positioning in urology and gynecology procedures with four-part leg section.



Advanced orthopedic and trauma applications supported by Extension Unit with MIS hip device.



Spine and neuro procedures on universal OR Table with modular structure and broad range of accessories.

Technical Details

	TruSystem™ 7000 U14	
Material-N°	TruSystem 7000 U14 (MBW) ¹ Pad TS7000 U14 B PadTS7000 U14 URO B	Mat.-No. 2065385 Mat.-No. 2065417 Mat.-No. 2065670
Height adjustment	59 cm – 114 cm / ± 1 cm (23" – 45" / ± 0.4")	
Trendelenburg/Anti-Trendelenburg	± 45°	
Tilt	± 30°	
Longitudinal slide	46 cm ± 1 cm (18" / ± 0.4")	
Joint leg section up/down	+ 90° / - 105°	
Joint back section up/down	± 90°	
Maximum load	450 kg (1000 lbs.)	

¹Motorized Base with Wireless Interface

Not all products / options are available in all countries. For further information about our products or services, please contact your local Trumpf Medical representative or visit our webpage:

www.trumpfmedical.com

Trumpf Medical reserves the right to make changes without notice in design, specifications and models. The only warranty Trumpf Medical makes is the express written warranty extended on the sale or rental of its products. ©2016 TRUMPF Medizin Systeme GmbH + Co. KG, ALL RIGHTS RESERVED.

TRUMPF Medizin Systeme GmbH + Co. KG
Carl-Zeiss-Straße 7-9 · 07318 Salfeld · Telefon +49 (0) 3671 586-0 · E-Mail info@trumpfmedical.com

Flyer ID: 2065954, 22-SEPT-2016